



Superior High Consequence Infectious Disease Protection

ViroGuard® 1 Designed to protect healthcare workers that are evaluating PUIs for Ebola or other HCIDs who are clinically stable and do not have bleeding, vomiting, or diarrhea, i.e. a **dry patient**.

ViroGuard® 1 fabric passes ASTM F1670 and F1671 for blood and blood-borne pathogens. Taped seams for additional protection.

Available Garments




Coverall with Boots Only. Elastic Wrists. Thumb Loops at Cuff. Zipper Front with Sealable Storm Flap. Taped Seams.

#W2401
Medium - 4X
25 per case



Coverall with Attached Hood & Boots. Elastic Wrists. Thumb Loops at Cuff. Zipper Front with Sealable Storm Flap. Taped Seams.

#W2404
Medium - 4X
25 per case



Coverall with Hood Only. Elastic Wrists. Thumb Loops at Cuff. Zipper Front with Sealable Storm Flap. Taped Seams.

#W2407
Medium - 4X
25 per case



Hood with Elastic Face, Tie Strings, Taped Seams.

#W2463
100 per case



Boot Cover with Elastic Top & Ankle. Skid Resistant Sole. Taped Seams.

#W2405-L
Height: 17"
50 per case

#W2406-XL
Height: 23"
50 per case

Available Garments



Coverall with Boots Only. Elastic Wrists. Thumb & Middle Finger Loops at Cuff. Zipper Front with Sealable Storm Flap. Double-Sided Tape at Wrists for Glove Adhesion. Mandarin Collar with Velcro Closure at Zipper. Taped Seams.

#W2501
Medium - 4X
25 per case



Coverall with Elastic Wrists & Ankles. Thumb & Middle Finger Loops at Cuff. Zipper Front with Sealable Storm Flap. Double-Sided Tape at Wrists for Glove Adhesion. Mandarin Collar with Velcro Closure at Neck. Red Fabric Pull Tab on Zipper. Taped Seams.

#W2502
Medium - 4X
25 per case



Hood with Elastic Face, Tie Strings, Taped Seams.

#W2563
100 per case



Boot Cover with Elastic Top & Ankle. Skid Resistant Sole. Taped Seams.

#W2505-L
Height: 17"
50 per case

#W2506-XL
Height: 23"
50 per case

ViroGuard® 2 Features Developed With Healthcare Workers' Input

Improves workers' safety and supports donning & doffing process.



Thumb & Finger Loops
Prevents sleeves from riding up.



Double Sided Tape Around the Wrists:
Helps keep gloves in place.



Large Zipper Pull with Red Fabric Pull Tab:
Easier to see and use!



Mandarin Collar with Velcro Closure:
Adjustable collar creates a better seal around the neck.



Zipper Front with Sealable Storm Flap has StarterTab for Tape:
Easier for a gloved hand to pull and expose the tape for sealing storm flap.



Visible Sizing Tag:
Helps materials management teams identify size easily once suits removed from box.

Agent Testing



AGENT	TEST METHOD	COMPONENT	RESULT
Blood	ASTM F1670	Fabric	Pass
Blood-Borne Pathogens	ASTM F1671	Fabric	Pass



AGENT	TEST METHOD	COMPONENT	RESULT
Blood	ASTM F1670	Fabric	Pass
Blood-Borne Pathogens	ASTM F1671	Fabric	Pass
Blood	ASTM F1670	Seams	Pass
Blood-Borne Pathogens	ASTM F1671	Seams	Pass

